



- SSMMS
- Melt Blown
- Active PP 25
- Melt Blown
- SSMMS

**Anti-Viral | Anti-Bacterial**

SaniShield is an **ACTIVE N95** facemask, consisting of 5 layer for Non-Woven fabric. The **ACTIVE N95** facemask is manufactured using Polyhexamethylene biguanide (PHMB) treated fabric which is proven to kill 99% of bacteria and virus (including Covid-19), in tests conducted by leading research labs globally.

SaniShield, an **ACTIVE N95** facemask possess **Anti-Viral & Anti-Bacterial** properties and has been tested as per ISO 20743:2013 | ISO 18184:2019 standards.

SaniShield an **ACTIVE N95** facemask, retains its activity for at least 30 cleaning cycles.

## SALIENT FEATURES

- Soft, Lightweight
- MULTIPURPOSE USE: Covid-19 Pandemic, Offices, Hospital, Pharma, Medical Stores, Laboratory, Food Sectors, Household, Industry and other purposes.

## 5 LAYERS OF FILTRATION

- SSMMS (Non-Woven 25 GSM) layer prevents entry of droplets into nose and mouth.
- Melt Blown (Non-Woven 25 GSM) layer filters bacteria, virus (including Covid-19 Virus), dust, allergens.
- PP (Non-Woven 25 GSM) layer treated with Anti-Bacterial and Anti-Viral solution.
- Melt Blown (Non-Woven 25 GSM) layer filters bacteria, virus (including Covid-19 Virus), dust, allergens
- SSMMS (Non-Woven 25 GSM) layer prevents entry of droplets into nose and mouth

## Anti-Viral | Anti-Bacterial



### TECHNICAL DETAILS:

- PRODUCT : SaniShield [ACTIVE N95 Face Mask]
- TESTED AS PER ISO 20743:2013 & ISO 18184:2019 STANDARD
- SITRA reports of all 5 fabric layers available

### INSTRUCTIONS

- Always store your **ACTIVE N95** mask in original packaging
- Wash your hand before & after handling your **ACTIVE N95** mask
- To clean dust and debris, use a moistened cloth on inner most and outer most layer of your **ACTIVE N95** mask.
- Replace your **ACTIVE N95** mask when it becomes dirty or gets damaged

### ⚠ WARNING

- Do not wear your mask while sleeping or doing strenuous exercise
- Do not share the mask with any other individual
- Discontinue to use your mask, in case you experience trouble breathing
- Do not dry clean or use washing machine or tumble dry for cleaning your mask



**Biosense Technologies Pvt. Ltd.**

A-233, Road No. 21-Y, Wagle Industrial Estate, Thane, Maharashtra- 400604  
Phone: 1800-123-7890 | Email: [contact@biosense.in](mailto:contact@biosense.in) | Website: [www.biosense.in](http://www.biosense.in)

**Biosense**  
A Tulip Diagnostics Company